

WATER

It is extremely important that appropriate water is used in sterilisers. Using tap or tank water will void your autoclaves warranty as well as leading to a build up of scale in the chamber, valves and pipe-work, which in turn can lead to unwanted breakdowns. The Australian Standard AS/NZS 4815-2006 requires that autoclave water reservoirs are drained weekly.

DISTILLERS



Produce your own autoclave quality distilled water with a **Cominox** water distiller.

Cominox Water Distiller, 230mm diameter, 360mm high, 4 litre capacity, 3.5kg, 580 Watts.

Part No. W-DIST.....\$490.00

REVERSE OSMOSIS



The **AquaClave RO P5** avoids the administrative inefficiency and high cost of purchasing bottled distilled water. This Australian developed process produces up to 50 litres per day of water with less than 1ppm in total dissolved solids. The AquaClave is supplied with a Total Dissolved Solids meter to check that the treated water meets the required high quality.

This system may be plumbed directly into the steam sterilisation equipment. This allows the water to be diverted to a separate faucet providing abundant water for the thorough rinsing of instruments prior to their placement into the autoclave.

AquaClave RO P5 Reverse Osmosis SystemPart No. W-ROP5\$1125.00

See our website www.autoclave.com.au for Australia's widest range of autoclaves. We have autoclaves to suit all needs including: beauty, day surgeries, dental, medical, laboratory, hospital, orthodontic, industrial, podiatry and veterinary.